
PM-07716


FemoSeal™
Vascular Closure System

INSTRUCTIONS FOR USE



1a


1b


2


3a


3b


4a


4b


4c

4d


5a


5b


6

7

8

9

ENGLISH

Device description


FemoSeal™ Vascular Closure System, manufactured by Terumo Medical Corporation (TMC), is a resorbable vascular closure device designed to achieve femoral arterial hemostasis after percutaneous catheterization through the common femoral artery.

FemoSeal closure elements consist of two resorbable polymer discs, the inner seal and the outer locking disc which are held together by a resorbable multifilament, thereby mechanically sealing the arteriotomy. After being deployed through the cone housing sheath, the outer locking disc is tamped onto the multifilament shaft resulting in mechanical closure of the puncture site between the inner and outer locking disc. The inner seal and outer locking disc are held in place via friction on the multifilament shaft. Hemostasis is achieved by mechanical means. The closure elements are degraded by hydrolysis.

The implantable closure components, inner seal and outer locking disc, are degraded in eighteen (18) months while the multifilament is estimated to be degraded and absorbed via hydrolysis by surrounding tissue after two (2) to three (3) years. The ability of the implantable closure components to maintain compression on the arteriotomy has been tested at up to eight (8) hours. The degradation products are metabolized and excreted in the urine or expired as carbon dioxide via the lungs. No accumulation effects have been observed in animal studies. The additional components of FemoSeal that will come in contact with living tissue or blood during the procedure are the cone housing sheath, dilator, tamping tube, pusher and the guidewire.

The FemoSeal unit is packaged together with a dilator and a 0.038" (0.97 mm) OD, 27.5" (70 cm) length guidewire and a J-straightener. The guidewire is contained within a polyethylene tube. All parts are packaged in a fixed paper tray. The packaging also includes the patient implant card. The FemoSeal Vascular Closure System does not require additional accessories for the device to function as intended and the device is only intended to be used with the dilator and guidewire included within the product's packaging.

FemoSeal™ VCS consists of the following components and materials.

Part Number	Description	Materials
1	0.038" (0.97 mm) Guidewire with a guidewire J-straightener	Guidewire: Stainless Steel J-Straightener: Polypropylene, white pigment Clip: Polyethylene Tubing: Polyethylene
2	FemoSeal Dilator	Hub: Tetrahydrofuran, Blue pigment, Polybutylene terephthalate (PBT) Tube: High density polyethylene (HDPE), black pigment Lubricant: Silicone oil
3/4 FemoSeal Unit	Molded RD7 – Inner Seal and Outer locking Disc *	Copolymer between glycolide, trimethylene carbonate, ε-caprolactone, and TMP: trimethylolpropane (initiator)
	Multifilament *	Segmented copolymer between L,L-Lactide, trimethylene carbonate, ε-caprolactone, and 1,3 propanediol (initiator) Coating is copolymer between glycolide, ε-caprolactone, and L-lysine
	Tamping Tube	Polypropylene
	Pusher	Stainless steel
	Slider	Polybutylene terephthalate (PBT)
	Housing	Polybutylene terephthalate (PBT)
	Button	Polybutylene terephthalate (PBT)
	Button Lid	Polyoxymethylene (POM)
	Sleeve	Polybutylene terephthalate (PBT)
	Spacer	Polybutylene terephthalate (PBT)
	Cone	Polypropylene
	Housing Lid	Polybutylene terephthalate (PBT)
	Safety Catch	Polybutylene terephthalate (PBT), blue pigment
	Sleeve Lid	Polybutylene terephthalate (PBT)
	Tube Gasket	Silicone shore A 70
	Housing Gasket	Silicone shore A 70
	Spring	Stainless Steel
	Cone Housing Sheath	Polypropylene Silicone Dow Corning 360, Hexane

* Implantable portion of device; Implantable components are MRI safe

FemoSeal™ VCS is sterilized via Ethylene Oxide (EO) and is a single-use medical device. The FemoSeal VCS device is designed in such a manner, that it cannot be re-used.

Mode of Action:

The FemoSeal™ VCS promotes mechanical hemostasis by deploying implantable closure elements inside and outside the common femoral artery. The mechanism that prevents bleeding is provided by the inner seal which seals the puncture site, and the outer locking disc which holds the inner seal in position.

Clinical Benefits

FemoSeal VCS provides the following clinical benefits relative to vascular closure after common femoral arterial puncture:

- Reduction in time to hemostasis¹
- Reduction in time to ambulation¹

1. Cox, T., Blair, L., Huntington, C., Lincourt, A., Sing, R., & Heniford, B. T. (2015). Systematic Review of Randomized Controlled Trials Comparing Manual Compression to Vascular Closure Devices for Diagnostic and Therapeutic Arterial Procedures. *Surgical technology international*, 27, 32-44.

Intended use

FemoSeal Vascular Closure System is a medical device intended for closure of an arterial puncture after percutaneous catheterization through the common femoral artery.

Intended users

FemoSeal Vascular Closure System's intended users are physicians with training qualifying them to perform arterial access and closure for endovascular procedures through the common femoral artery and have participated in a Terumo Medical Corporation FemoSeal physician instruction program.

Indication for use

FemoSeal Vascular Closure System is indicated for use in closing the common femoral arterial puncture (arteriotomy) in patients who have undergone percutaneous catheterization using a 7F (2.33 mm) or smaller procedural sheath.

Contraindications

FemoSeal Vascular Closure System is contraindicated in patients with arteriotomies in which sheaths or devices larger than 7F (2.33 mm) have been used.

Target patient populations

The safety and effectiveness of FemoSeal Vascular Closure System has been established in the patients 18 years of age and older who have undergone percutaneous catheterization using a 7F (2.33 mm) or smaller procedural sheath.

Special patient populations

The safety and effectiveness of FemoSeal Vascular Closure System has not been established in the following patient populations:

- Patients with pre-existing autoimmune disease.
- Patients undergoing therapeutic thrombolysis.
- Patients with clinically significant peripheral vascular disease.
- Patients with uncontrolled hypertension (> 220 mmHg systolic or > 110 mmHg diastolic).
- Patients with a bleeding disorder, including thrombocytopenia (< 100,000 platelet count), or anemia (Hgb < 10 mg/dl).
- Patients having an inner lumen of the common femoral artery smaller than 5 mm.
- Patients with a vascular graft or stent at the puncture site.
- Patients who are pregnant or lactating.
- Patients below the age of 18 years.
- Patients with renal and/or hepatic impairment.
- Patients with other relevant co-morbidity.
- Population with specific racial and/or ethnic origins.
- Patients with myocardial infarction within 72 hours.

Warnings

- Do not use if the package has been damaged or any sterile barrier is not intact.
- Do not use after expiry date – the biodegradable components may not perform adequately.
- Do not use if any items appear damaged or defective in anyway.
- Use of FemoSeal Vascular Closure System where bacterial contamination of the procedural sheath or the surrounding tissue may have occurred, may cause infection.
- If it is suspected that the posterior arterial wall has been punctured or more than one arterial puncture has been made, do not rely solely on FemoSeal Vascular Closure System to achieve arterial hemostasis. Use additional manual or mechanical compression.
- If the puncture site is at, or distal to the bifurcation of the femoral artery, FemoSeal Vascular Closure System should not be used due to the risk of the Inner Seal being positioned incorrectly. This event may result in bleeding complications and/or disruption to normal blood flow.
- If there is persistent arterial bleeding from the incision site, do not cut the multifilament until hemostasis is achieved. If hemostasis is not achieved, fasten the multifilament with a sterile wound dressing and apply supplementary compression until hemostasis is achieved. In the case of persistent arterial bleeding, significant bleeding complications may occur which could result in patient injury or death.
- Do not use the FemoSeal Vascular Closure System if the puncture site is proximal to the inguinal ligament as this may result in a retroperitoneal hematoma.
- Patients at a higher risk for bleeding may suffer increased blood loss, requiring a transfusion.
- If the inner disc becomes detached it can potentially cause a thrombotic and/or embolic event.

Precautions

- FemoSeal Vascular Closure System deployment procedure should be performed by physicians/clinicians with adequate training in the use of the device.
- Perform a limited femoral angiogram or an ultrasound-guided femoral puncture prior to FemoSeal Vascular Closure System deployment.
- Discontinue procedure if:
 - Lumen diameter of common femoral artery < 5 mm.
 - Stenosis and/or significant plaque present in the vicinity of the femoral arterial puncture site.
 - Arterial puncture is at, or distal to, the common femoral artery bifurcation.
 - Anomalous branches or vessel abnormalities present in the vicinity of the femoral arterial puncture site.
- If any hematoma is present, extra care must be taken for correct insertion of the cone housing sheath into the artery.
- Use a single wall puncture technique. Do not puncture the posterior wall of the artery.
- Observe sterile technique at all times when using FemoSeal Vascular Closure System.
- FemoSeal Vascular Closure System is for single use only and should not be resterilized or reused in any manner. The FemoSeal unit is designed in such a manner, that it cannot be re-used.
- If deployment of the Inner Seal meets unexpected resistance, discontinue the procedure.
- For correct deployment of Inner Seal, reposition your thumb so that the button can spring back freely.
- Ensure the tip of the cone housing sheath of FemoSeal Unit is under the skin surface before deployment of the Outer Locking Disc. The Outer Locking Disc may inadvertently be deployed above skin level in patients with a short distance between the common femoral artery and the skin level.
- For correct deployment, a skin incision may be necessary before deployment.
- If the Inner Seal is inside the artery, but it is not possible to fully depress the button and deploy the Outer Locking Disc, surgery may be required to remove FemoSeal Closure Elements.
- If repuncture of the same femoral artery becomes necessary within 18 months, repuncture should be made at least one centimeter proximal to the previous FemoSeal Vascular Closure System access site.
- Instruct the patient to follow physician's orders regarding closure site inspection.
- Instruct the patient to carry the Patient Information Card for the next 18 months.

Precautions at time of discharge

Before considering discharge, assess the patient for the following clinical conditions:

- Bleeding and/or hematoma at the closure site.
- Pain while walking.
- Signs of infection at the closure site.

Adverse Events, Harms, and Complications

Known or foreseeable adverse events, harms, and complications associated with the use of FemoSeal VCS include:

- Allergic Reaction
- Aneurysm
- AV Fistula
- Blood Loss/Bleeding
- Death
- Device Failure
- Ecchymosis
- Embolism
- Foreign Body Reaction
- Hematoma
- Hemorrhage
- Infection
- Inflammatory Reaction
- Numbness
- Pain
- Patient Discomfort
- Procedure delay
- Pseudoaneurysm
- Retroperitoneal Bleed
- Sepsis
- Thromboembolism
- Thrombosis
- Vessel occlusion/Lower Limb ischemia
- Vessel Perforation
- Vessel Tissue Dissection/Laceration

Directions for Use

The procedure is described in the following steps:

1. Assess puncture site location by fluoroscopy or ultrasound.

Perform a limited femoral angiogram or an ultrasound-guided femoral puncture prior to FemoSeal Vascular Closure System deployment.

CAUTION: Discontinue procedure if:

- Lumen diameter of common femoral artery < 5 mm.
- Stenosis and/or significant plaque present in the vicinity of the femoral arterial puncture site.

- Arterial puncture is at, or distal to, the common femoral artery bifurcation.

- Anomalous branches or vessel abnormalities present in the vicinity of the femoral arterial puncture site.

CAUTION: If any hematoma is present, extra care must be taken for correct

insertion of the cone housing sheath into the artery.

2. Position the cone housing sheath of FemoSeal Unit inside the artery.

- Insert FemoSeal Dilator into FemoSeal Unit until they snap together and marked signs meet (see figures 1a, 1b).
- **NOTE:** Grip FemoSeal Dilator near its tip, advance carefully in small increments.
- Insert the guidewire through the procedural sheath.
- Remove the procedural sheath and hold the guidewire in position (see figure 2).
- Thread FemoSeal Unit over the guidewire and insert the cone housing sheath into the artery (see figures 3a, 3b).

NOTE: Leave approximately 25 cm of the guidewire visible.

3. Grip and squeeze the wings of FemoSeal Safety Catch, lift and retract FemoSeal Safety Catch slowly along with the guidewire until completely removed (see figures 4a, 4b, 4c, and 4d).

NOTE: Hold the device in position with the opposite hand. When FemoSeal Safety Catch with dilator has been retracted a few centimeters (2 - 4 cm), blood should appear in the proximal part of the cone housing sheath, verifying that the cone housing sheath is in the artery.

NOTE: If no blood appears, retraction may be disrupted. Reintroduce the dilator and then the cone housing sheath into the artery and repeat the steps from 3.

4. Deploy the Inner Seal by pressing down the button completely in one continuous movement (see figures 5a, 5b).

CAUTION: If deployment of the Inner Seal meets unexpected resistance, discontinue the procedure.

5. Pull back FemoSeal Unit until the depressed button springs back (see figure 6).

CAUTION: For correct deployment of Inner Seal, reposition your thumb so that the button can spring back freely.

NOTE: The button springs back once the Inner Seal makes contact with the inside of the arterial wall and a resistance is generated.

NOTE: The Inner Seal brings about the hemostasis.

6. Deploy the Outer Locking Disc by completely pressing down the button once again (see figure 7).

CAUTION: Ensure that the tip of the cone housing sheath of FemoSeal Unit is under the skin surface before deployment of the Outer Locking Disc. The Outer Locking Disc may inadvertently be deployed above skin level in patients with a short distance between the femoral artery and the skin level. For correct deployment, a skin incision may be necessary before deployment.

NOTE: A steep puncture angle shortens the distance between the femoral artery and the skin level.

NOTE: Keep a slight tension in the multifilament during deployment of the Outer Locking Disc until the button is completely pressed down. This step locks the Inner Seal against the arterial wall.

CAUTION: If the Inner Seal is inside the artery, but it is not possible to fully depress the button and deploy the Outer Locking Disc, surgery may be required to remove FemoSeal Closure Elements.

7. Remove FemoSeal Unit (see figure 8).

8. Verify hemostasis.

WARNING: If there is persistent arterial bleeding from the incision site, do not cut the multifilament until hemostasis is achieved. If hemostasis is not achieved, fasten the multifilament with a sterile wound dressing and apply supplementary compression until hemostasis is achieved. In the case of persistent arterial bleeding, significant bleeding complications may occur which could result in patient injury or death.

9. Cut the multifilament below skin level using a sterile instrument (see figure 9).

Storage and packaging

- Sterile and non-pyrogenic in unopened and undamaged package.
- FemoSeal Vascular Closure System should be stored in a cold location (room temperature not exceeding 25 °C).
- FemoSeal Vascular Closure System contains resorbable materials that de-grade by exposure to heat and moisture; therefore, the device may not be re-sterilized.
- This device is sterilized by Ethylene Oxide (EO) and is provided in a sterile state and is a single-use medical device. The FemoSeal VCS device is designed in such a manner, that it cannot be re-used.
- After the packaging has been peeled open, verify the sterile barrier is intact by inspecting the packaging to ensure there are no channels, wrinkles, or voids in the seal of the packaging.

Disposal of the FemoSeal Vascular Closure System

This device should be disposed of using standard hospital procedures and universal precautions for biohazardous waste. All relevant laws and regulations governing such disposal should be followed.

Summary of Safety and Clinical Performance

The summary of safety and clinical performance (SSCP) is intended to enhance transparency by providing public access to an updated summary of clinical data and other information about the safety and clinical performance of the medical device. The SSCP is available in two versions, a version for healthcare professionals and a separate version for patients. The SSCP for the FemoSeal VCS medical device can be found at <https://ec.europa.eu/tools/eudamed>.

Table: FemoSeal Vascular Closure System Troubleshooting

Issues	Scenario	Possible Causes
Non-Deployment	The closure elements are pulled out from the artery, upon withdrawal of the device,	The inner seal was deployed outside the artery.
Unexpected Resistance	The inner seal meets unexpected resistance during deployment (pressing plunger)	The tip of the cone housing sheath may be impinging on the posterior wall of the artery.
The Plunger is blocked	The plunger does not spring back when pulling back	The thumb is on the plunger, blocking the plunger's correct ejection. FemoSeal VCS not in tension: the inner seal is not opposed to the vessel wall. Not enough space in the artery for the inner seal placement against the vessel wall. The device is damaged or defective..
Multifilament breakage	Upon withdrawal of the FemoSeal VCS, a multifilament breakage is suspected.	The device was not used according to the instructions for use. The device is damaged or defective.
Disc deployed above skin level	The outer locking disc is deployed above skin level	The outer locking disc may inadvertently be deployed above skin level in patients with a short distance between the femoral artery and the skin level. The inner seal was deployed outside the artery.
Persistent Arterial Bleeding	Persistent arterial bleeding from the incision site after deployment.	The posterior arterial wall has been punctured. More than one arterial puncture was made. The inner seal was deployed outside the artery.

Explanation of symbols:

Vascular Closure System	Vascular Closure System.
	Consult Instructions for use.
	Complies with the Medical Device Regulation 2017/745.
	Quantity.
	Store in a cool location (room temperature).
	Catalogue number.
	Use-by date.
	Sterilized using Ethylene Oxide.
	Fragile, handle with care.
	Batch code.
	Do not use if package is damaged.
	Do not resterilize.
	For single use only. Do not reuse.
	Keep dry.
	Keep away from sunlight.
	Manufacturer.
	Date of manufacture.
	Country of Origin
Made in USA	Instructions for Use..
INSTRUCTIONS FOR USE	

Possible Actions	Preventive Actions
Apply manual or mechanical pressure per standard procedure. Examine the device to ensure all bioabsorbable components have been withdrawn. Otherwise, surgery may be required to remove the closure elements.	Verify that the cone housing sheath is in the artery before deploying the inner seal (blood in chamber). Pull up the carrier to avoid a small angle versus skin ($\geq 45^\circ$)
FemoSeal VCS should not be used. Apply manual or mechanical pressure per standard procedure	If the puncture site is at, or distal to the bifurcation of the femoral artery, FemoSeal VCS should not be used due to the risk of the inner seal being positioned incorrectly.
-Make sure your thumb has been removed from the plunger, so it can spring back freely. Create a greater angulation of the FemoSeal ($\geq 45^\circ$) compared to skin level. The button springs back once the inner seal makes contact with the inside of the arterial wall and resistance generated.	No preventive action. Keep maintaining the FemoSeal VCS in tension with a greater angulation ($\geq 45^\circ$). Perform a femoral angiogram or ultrasound to ensure lumen diameter of the artery $\geq 5\text{mm}$ prior to use of FemoSeal VCS
Apply manual or mechanical pressure per standard procedure. Surgery may be required to remove the closure elements.	No preventive action
Do not cut the multifilament. Fasten the multifilament with a sterile wound dressing and apply supplementary compression. Surgery may be required.	Do not force, but only have tactile feedback of the inner seal against vessel wall. Ensure the tip of the FemoSeal cone housing sheath is under the skin surface before deployment of the outer locking disc. Verify that the cone housing sheath is in the artery before deploying the inner seal.
If there is persistent arterial bleeding from the incision site, do not cut the multifilament until hemostasis is achieved. If the hemostasis is not achieved, fasten the multifilament with a sterile wound dressing and apply supplementary compression until hemostasis is achieved. If it is suspected that the posterior arterial wall has been punctured or more than one arterial puncture has been made, do not use FemoSeal to achieve arterial hemostasis. Use additional manual pressure or a compression assist device. Verify that the cone housing sheath is in the artery before deploying the inner seal.	Perform a femoral angiogram or ultrasound to ensure lumen diameter of the vessel $\geq 5\text{mm}$ prior to use of FemoSeal VCS. Verify the access site closure via an ultrasound guidance.

Authorized Representative in the European Community.

EC REP

Medical device.

MD

Single sterile barrier system with protective packaging outside.

Single sterile barrier system with protecting packaging inside.

UDI

Unique Device Identification.

MR

MR - Implantable components are MRI safe

?

Patient name/identification.

+

Healthcare institution.

31

Implant date.

?

Patient information website.

?

Importer.

Guide Wire Length

Guide Wire Length

Max Guide Wire O.D.

Max Guide Wire Outer Diameter

Unique Device Identification - Device Identification.

UDI-DI

Terumo Europe N.V.

Interleuvenlaan 40

3001 Leuven Belgium

+32 16 38 12 11

Australian sponsor:

Terumo Australia Pty Ltd

Macquarie Park NSW 2113

Australia

Legal Manufacturer:

Terumo Medical Corporation

265 Davidson Avenue, Suite 320

Somerset, NJ 08873 USA

+1 732 302 4900

+1 800 283 7866

www.terumomedical.com/ifu

In Case of Adverse Events, Report to:

Report adverse events to the legal manufacturer, EC representative, and your EU Competent Authority of your member state in which the user and or patient is established.

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Unless otherwise noted,™ indicates a registered or unregistered trademark or service mark owned by, or licensed to, Terumo Europe N.V. or Terumo Medical Corporation. FemoSeal is a registered and/or unregistered trademark of Terumo Europe N.V. © 2018 Terumo Medical Corporation. All rights reserved.