Summary of Safety and Clinical Performance

FemoSeal™ Vascular Closure System

Manufacturer	Terumo Medical Corporation 265 Davidson Ave, Suite 320 Somerset, NJ 08873 USA
Manufacturer's SRN	US-MF-000019594
Basic UDI-DI	38970FS79
EMDN number	C900199 – Haemostasis Systems-Other
Device Class	Class III
Year Device was First CE Marked	2005
EU Representative	Terumo Europe, N.V. Interleuvenlaan 40, 3001 Leuven Belgium SRN: BE-AR-000001433
Notified Body	NSAI CE 0050

Table of Contents

1.0	Summary for Patients.....	4
1.1	Intended Use of the Device.....	4
1.1.1	Intended Purpose.....	4
1.1.2	Indications for Use.....	4
1.1.3	Contraindications	4
1.1.4	Intended patient population.....	4
1.2	Device Description.....	5
1.2.1	Description of the Device	5
1.2.2	Raw Materials & Components.....	6
1.3	Risks and Warnings	8
1.3.1	Remaining Risks and Undesirable Effects.....	8
1.3.2	Warnings and Precautions.....	9
1.3.3	Field Safety Corrective Actions.....	9
1.4	Summary of Clinical Evaluation and Post-Market Clinical Follow-Up.....	9
1.4.1	Clinical Background of the Device.....	9
1.4.2	Clinical Evidence Supporting CE Marking	10
1.4.3	Safety	20
1.5	Possible Therapeutic Alternatives.....	20
2.0	Bibliography	21
	Clinical Background Publications.....	21
	Publications Providing Safety and Performance Data on FemoSeal™ VCS	21
3.0	Revision history	23

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the FemoSeal™ Vascular Closure System (VCS).

A summary of the safety and clinical performance of the device, intended for patients, is provided in **Section 1.0**.

1.0 SUMMARY FOR PATIENTS

The information presented below is intended for patients or lay persons.

This SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation.

This SSCP is not intended to replace an Implant Card or the Instructions for Use to provide information on the safe use of the device.

1.1 Intended Use of the Device

1.1.1 Intended Purpose

FemoSeal™ Vascular Closure System is a medical device used to close a hole in an artery (common femoral arterial puncture) of your groin.

1.1.2 Indications for Use

FemoSeal™ Vascular Closure System is approved for use in closing a hole in the artery in the groin.

1.1.3 Contraindications

FemoSeal™ is not to be used (contraindicated) in patients with holes in the artery that are larger than 2.33 mm (7 French) diameter.

1.1.4 Intended patient population

FemoSeal™ is approved to be used in patients 18 years of age and older.

The use of FemoSeal™ has not been tested in the following patients:

- Patients getting blood thinning medication (therapeutic thrombolysis).
- Patients with narrow arteries (clinically significant peripheral vascular disease) at the puncture site.
- Patients with high blood pressure (Uncontrolled hypertension of > 220 mmHg first number in a blood pressure reading [systolic] or > 110 mmHg second number in a blood pressure reading [diastolic]).
- Patients with a bleeding disorder (thrombocytopenia) or low iron (anemia).
- Patients with small arteries (common femoral artery smaller than 5mm).
- Patients with kidney or liver disease (renal and/or hepatic impairment).


- Patients with other diseases (relevant co-morbidity).
- Population of specific race.
- Patients who had a heart attack (myocardial infarction) within 72 hours.
- Patients with implanted device at the hole in the artery (Stent or vascular graft).
- Patients who are having a baby or breast feeding.
- Patients below the age of 18 years.

1.2 Device Description

1.2.1 Description of the Device

FemoSeal™ is a device used to close a hole in the femoral artery (arterial puncture). FemoSeal™ implantable parts include two degradable plastic (polymer) discs that are left in the body and a suture. The discs make a sandwich like seal held together by suture (multifilament). See **Figure 1.2**.

Figure 1.1 FemoSeal™ VCS Unit with the Guidewire and Dilator

1. Guidewire 0.038" (0.97 mm)

2.Artery stretching tool (Dilator)


3.FemoSeal™ Unit

4.Tube that allows passage in and out of artery (Sheath)

Figure 1.2 Inner Seal and Outer Locking Disc, Respectively

Figure 1.3 Inner Seal and Outer Locking Disc Deployed

The discs left in the body dissolve within eighteen (18) months. The suture (multifilament) dissolves within two (2) to three (3) years.

1.2.2 Raw Materials & Components

FemoSeal™ parts that are left in the body will dissolve (degrade) over time. Pictures in **Figure 1.3**. **Table 1.1** below lists the materials in the device.

Table 1.1 Raw Materials in FemoSeal™ VCS

Part Number	Description	Materials
1	0.038" (0.97 mm) Guidewire with a guidewire J-straightener	Guidewire: Stainless Steel J-Straightener: Polypropylene, white pigment Clip: Polyethylene Tubing: Polyethylene
2	FemoSeal™ Dilator	Hub: Tetrahydrofuran, Blue pigment, Polybutylene terephthalate (PBT)

			Tube: High density polyethylene (HDPE), black pigment Lubricant: Silicone oil
3/4	FemoSeal™ Unit	Molded RD7 – Inner Seal and Outer locking Disc *	Copolymer between glycolide, trimethylene carbonate, ε-caprolactone, and TMP: trimethylolpropane (initiator)
		Multifilament *	Segmented copolymer between L,L-Lactide, trimethylene carbonate, ε-caprolactone, and 1,3 propanediol (initiator) Coating is copolymer between glycolide, ε-caprolactone, and L-lysine
		Tamping Tube	Polypropylene
		Pusher	Stainless steel
		Slider	Polybutylene terephthalate (PBT)
		Housing	Polybutylene terephthalate (PBT)
		Button	Polybutylene terephthalate (PBT)
		Button Lid	Polyoxymethylene (POM)
		Sleeve	Polybutylene terephthalate (PBT),
		Spacer	Polybutylene terephthalate (PBT)
		Cone	Polypropylene
		Housing Lid	Polybutylene terephthalate (PBT)
		Safety Catch	Polybutylene terephthalate (PBT), blue pigment
		Sleeve Lid	Polybutylene terephthalate (PBT)
		Tube Gasket	Silicone shore A 70
		Housing Gasket	Silicone shore A 70
		Spring	Stainless Steel
		Cone Housing Sheath	Polypropylene Silicone Dow Corning 360, Hexane

* Implantable portion of device; implantable components are MRI safe

1.3 Risks and Warnings

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed. Your doctor has considered the risks and decided that this device is appropriate for your procedure.

1.3.1 Remaining Risks and Undesirable Effects

Harms associated with the use of the FemoSeal™ VCS as a result of residual risks are listed in **Table 1.2**. The occurrence rate is based on data collected from January 2015 through December 2021.

Table 1.2 Harms Associated with Residual Risks

Known or Foreseeable Harm	Occurrence Rate (%)	Mitigating Factors
Allergic reaction	0	Actions and plans to reduce risk as low as possible have been addressed by product design features, manufacturing guidelines, product labeling, and physician training.
Weakening of the artery wall (Aneurysm)	0	
Abnormal connection between an artery and a vein (AV Fistula)	7.68036E-05	
Blood loss / bleeding	0.05767947	
Death	0	
Bleeding under the skin (Ecchymosis)	0	
Blockage of an artery (Embolism)	0.000153607	
Irritation due to the device (Foreign body reaction)	7.68036E-05	
Blood Loss (Hematoma)	0.00430100	
Release of blood from a broken blood vessel (Hemorrhage)	0.00023041	
Infection	0.000230411	
Irritation, reddened area (Inflammation)	0	
Loss of sensation (Numbness)	0	
Pain	0.000230411	
Patient discomfort	0	

Known or Foreseeable Harm	Occurrence Rate (%)	Mitigating Factors
Procedure delay	0.018509658	
An incomplete healing or sealing of the artery where the hole was (Pseudoaneurysm/false aneurysm)	0.00384018	
Bleeding into the abdominal space (Retroperitoneal Bleed)	0.000460821	
Infection (Sepsis)	0	
Blockage of an artery due to a clot (Thromboembolism)	0.000153607	
Blockage of an artery due to a clot (Thrombosis)	0	
Blockage of an artery (Vessel occlusion or lower limb ischemia)	0.00921643	
Tear in the artery (Vessel perforation)	0	
Tear in the artery (Vessel tissue dissection/laceration)	0	

1.3.2 Warnings and Precautions

- Follow doctor's orders.
- Call your doctor if you have bleeding, pain, or signs of infection. Signs of infection are redness, fever, discharge and/or warmth at the site. Signs of a hematoma are swelling, lump in groin.
- Carry the Patient Information Card for the next 18 months.

1.3.3 Field Safety Corrective Actions

There was one recall involving FemoSeal™. The recall was due to an increase in product complaints. The failure was caused during manufacturing and has been corrected.

1.4 Summary of Clinical Evaluation and Post-Market Clinical Follow-Up

1.4.1 Clinical Background of the Device

FemoSeal™ has been on the market for 15 years.

1.4.2 Clinical Evidence Supporting CE Marking

FemoSeal™ approval is based on a study of 80 patients. The study showed that FemoSeal™ VCS was able to safely close the artery in at least 95% of patients.

Table 1.3 is a list of studies with FemoSeal™ VCS. The original version of the device was on the market until 2006. The current device version is in the marketplace today.

Table 1.3 Summary of Pivotal Studies

Pre-Market Functional Animal Testing	Device Version
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
First In Human Clinical Study	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (original device)
Completed Post-Market Clinical Studies	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (current device)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (current device)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (current device)
Planned / Ongoing Post-Market Clinical Studies	
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe	C11202 (current device)
Post Market Clinical Follow Up Activities	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (current device)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (current device)

Table 1.4 below describes all studies for FemoSeal™ VCS.

Table 1.4 Summaries of Pre-market Studies

Pre-Market Functional Animal Testing					
Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	2000 Radi Medical System AB Performed by Biomatech, France	Animal study; Exploratory investigation	The aim of the study was to evaluate the biocompatibility and biofunctionality of a biodegradable implant designed to achieve hemostasis in the vascular system. Report addresses reserve animals from Study Number 862 (Document # 1051-01)	N = 3 adult female sheep, 1 animal per time point Follow up for histopathologic analysis, macroscopic grading, and ultrastructural analysis • 12 months • 15 months • 18 months Color doppler ultrasound follow up – 30 days at carotid and femoral arterial access sites	Arterial patency • Post deployment blood flow assessment – 2/3 blood flow reduction associated with temporary arterial vasospasm that was attributed to procedure, both occurred at left common carotid artery • Color doppler ultrasound assessment (4 weeks) – 3/3 (100%) without major blood flow reduction 30 days after implantation Macroscopic observations • 12 months – patency confirmed macroscopically by absence of anatomically visible occlusion; no macroscopic local intolerance lesions (hemorrhage, necrosis, or neovascularization); encapsulation of internal disc appeared to be marked • 15 months - macroscopic observation showed no signs of inflammation; no residues of inner or outer discs identified; encapsulation tissue located in area of outer disc appeared the same as 12-month observation • 18 months – macroscopic observation showed no signs of local intolerance, outer disc visible (3/3 cases), degradation of inner disc complete; 1/3 (Left common carotid artery) – depression of artery wall observed at implant site (could be related to device deployment) Conclusions • Macroscopically, 12, 15 and 18 months after

					implantation, no necrotic, degenerative, or thrombotic signs were detected. - The degradation process of the inner and outer discs appeared to progress between 12 and 18 months and seemed complete at 18 months, with a whitish tissue covering the implanted area. - Histologically, after 18 months, polymer material was no longer visible but an endoluminal fibrous tissue was present. No necrotic, degenerative, or thrombotic lesions were noted.
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	2002 Radi Medical System AB Performed by Biomatech, France	Animal study; Exploratory investigation	The aim of the study was to evaluate the biocompatibility of a biodegradable implant designed to achieve hemostasis in the vascular system.	Follow up (after implantation) 12 weeks 24 weeks 36 weeks N = 12 sheep, 3 animals per time point, 3 reserve animals sacrificed at 12, 15, and 18 months Implantation sites – right and left femoral arteries, right and left common carotid arteries Color doppler ultrasound analysis	<u>Patency of carotid and femoral arteries after 4 weeks</u> - Major blood flow reduction 0/12 (0%) <u>No major blood flow reduction</u> - Severe blood flow impairment related to device insertion 0/12 (0%) - Blood flow reduction in left common carotid 2/12 (16.6%) due to temporary arterial vasospasm (attributed to surgery, not device) <u>Blood flow impairment by follow up</u> 12 weeks – 0/12 (0%) 24 weeks – 0/12 (0%) 36 weeks – 0/12 (0%) Blood pressure measurements (3 sheep – 36 weeks) - No major blood flow reduction detected Macroscopic data for integration and degradation (Grade 0-4) <u>Carotid sites</u>

					<i>Internal side</i> -encapsulation 12 week 1.3 24 week 0.8 36 week 0.7 -material degradation 12 week 2.2 24 week 2.3 36 week 0.7 <i>External side</i> -encapsulation 12 week 2.0 24 week 1.0 36 week 0.8 -material degradation 12 week 2.0 24 week 2.7 36 week 3.2 <u><i>Femoral sites</i></u> <i>External side</i> -encapsulation 12 week 1.0 24 week 2.7 36 week 1.0 -material degradation 12 week 2.0 24 week 2.4 36 week 3.0 <i>Internal side</i> -encapsulation 12 week 2.0 24 week 1.2 36 week 1.0 -material degradation 12 week 1.0 24 week 2.2 36 week 3.0 <u>Histological data for integration, degradation, and inflammatory parameters (mean of relevant values, Index 0-4) n=6</u> <u><i>Carotid sites</i></u> External fibrosis 12 weeks – 2.3 24 weeks – 2.3 36 weeks – 1.8

					Internal neointima 12 weeks – 3.0 24 weeks – 2.7 36 weeks – 3.3 Material degradation 12 weeks – 1.7 24 weeks – 3.0 36 weeks – 3 + Inflammatory parameters (macrophages) 12 weeks – 3.0 24 weeks – 1.7 36 weeks – 1.7 <u>Femoral sites</u> External fibrosis 12 weeks – 1.3 24 weeks – 1.5 36 weeks – 1.7 Internal neointima 12 weeks – 3.7 24 weeks – 3.3 36 weeks – 3.8 Material degradation 12 weeks – 2.0 24 weeks – 3.0 36 weeks – 3+ Inflammatory parameters (macrophages) 12 weeks – 3.5 24 weeks – 1.8 36 weeks – 1.8 <u>Histomorphological patency values</u> Carotid sites 12 week 74.0% 24 week 81.9% 36 weeks -89.43% Femoral sites 12 week – 42.2% 24 week – 73.21% 36 weeks – 78.71% <u>US vessel patency (mean %)</u> Carotid site 12 week 94.6% 24 week – 96.7% 36 week – 94.7% Femoral sites 12 week – 66.9%

					24 week – 63.6% 36 week – 75.1% Conclusions - No significant local intolerance sign was detected (absence of visible inflammatory, necrotic, or degenerative lesions) at 12, 24, or 36 weeks - No significant adverse tissular reaction recorded with regard to host and device - Degradation of implant was nearly complete at 36 weeks
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	2002 Radi Medical System AB Performed by Biomatech, France	Animal study	The purpose of the study was to evaluate the biocompatibility and biofunctionality of a biodegradable implant designed to achieve haemostasis in the vascular system. The tested device was intended to enter the puncture hole in arteries to achieve haemostasis. One disc entered the puncture hole in the artery and the other tightened the internal disc from outside the artery. Both discs were inserted through an introducer sheath. The discs were designed to prevent bleeding, to heal into the vessel wall and finally degrade and be absorbed. The study was designed to collect and analyze data after functional implantation in the sheep and pig to evaluate the capacity of the implant to maintain physiologic function	Implantation period / follow up is 5 weeks for all animals 4 sheep, 1 pig 14 closure elements 8 closure elements implanted in 2 sheep, 2 closure elements implanted in 1 pig When two closure elements were implanted into same femoral artery, they were implanted at a distance > 30 mm	Any adverse signs criteria (inflammation, necrosis, hemorrhage, or any other lesion) or any thin pseudo-intimal layer formation were recorded. Particular attention was devoted to initial tissular integration of the device, with a special emphasis on the internal disc encapsulation. Macroscopically, no significant signs of local intolerance or thrombus formation were observed. The macroscopic observation did not reveal any difference between the sheep or the pig after implantation of the closure elements in the femoral arteries. Macroscopically, no significant signs of local intolerance or thrombus formation were observed. The macroscopic observation did not reveal any difference between the sheep or the pig after implantation of the closure

			when used in the circulatory system, to determine the response of the host and the response of the device (physical integrity, tissular encapsulation at sacrifice) and to evaluate the degradation of the implant		elements in the femoral arteries. The histopathological analysis of the implantation site in sheep showed no local adverse reactions. The inner and outer discs showed signs of degradation and were in some cases broken into smaller parts. The inner discs were completely integrated into neointimal tissue. In the pig the histopathological analysis showed slightly more inflammatory signs around the inner and outer discs. In both animals, the multifilament exhibited slightly more inflammatory signs than in the disc material, which is normal for braided sutures due to the larger surface area. The vessel patency was measured after 5 weeks using three different methods, histomorphometric, histomorphometric in combination with measurement of vessel diameter by use of caliper and by the use of color doppler ultrasonography. The mean vessel patency in sheep was between 60 and 80 % depending on method (12 observations). In the pig, the mean vessel patency was between 40 and 80 % (2 observations).
First In Human Clinical Study					

Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	2004 Radi Medical Systems AB	Open label, single arm, multi-centre trial	The objectives of the study were to determine the safety, performance, and effectiveness of the FemoSeal Vascular Closure System (FemoSeal System) in general, and specifically the Closure Elements, in patients undergoing cardiac diagnostic or interventional cardiology procedures using femoral artery access.	2-part study Total n = 80 Part 1 n = 20 Part 2 n = 60 Indications: Part 1 – diagnostic cardiac catheterization using femoral artery as site of access Part 2 – patients undergoing diagnostic and interventional cardiac procedures using femoral artery as the site of access 30 day follow up	Outcomes <u>-hemostasis achieved with FS without compression</u> - Part 1 20/20 (100%) - Part2 58/60 (97%) Median time to hemostasis at time of procedure (part 2) – 1.0 minutes (Range immediate – 8 minutes) <u>Time to ambulation (minutes)</u> Part 1 – Median: 33.0 mean: 67.8 Part 2 – median: 45 mean: 87.5 <u>Incidence of major vascular complications – 30 day follow up</u> Part 1 -Free of vascular complications 20/20 (100%) Part 2 Pseudoaneurysm - 1/59 (1.69%) <u>Treatment failure</u> Part 1 – 0/20 (0%) Part 2 – 3/60 (5%) • 2 due to the failure to achieve hemostasis without use of compression • 1 due to major vascular complication, pseudoaneurysm <u>Persistent oozing</u>

					Part 1 – 0/20 (0%) Part 2- 8/60 (13%) <u>Hematoma</u> Part 1 – 0/20 (0%) Part 2 – 6/60 (10%)

Table 1.5 below summarizes all studies completed once the device went to market.

Table 1.5 Post-Market Studies

Completed Post-Market Clinical Studies					
Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
PMS FemoSeal® VCS #11202 (R1756-01) (2006) User Acceptability of New Packaging (SIS)	2006 RADI	Post Market Study, User Acceptability of New Packaging	Post-market study carried out on the model update FemoSeal Vascular® Closure System #11202 which has a modified delivery system compared to the first model, #11200, to document functionality. The objective of the study was to get feedback on the revised Instructions for Use (IFU), the training project for physicians and nurse, to confirm the assessments made in the risk analysis and to give input to any product quality improvements.	100 patients 63 male / 37 females - Mean age 68 years for female, 65 years for males 2 Swedish hospitals, 6 physicians 112 devices Data captured on patient record forms (CRF) Planned observation time – 30 minutes from time of application	Outcomes 1. Illustrations/Instructions in IFU considered informative and adequate 2. Immediate hemostasis (91/100)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	Dr. Nicolas Moes Innsbruck Medical University, Austria	Post -Market Confirmatory Assessment	The objective of this study was to assess safety and performance of FemoSeal® used during normal clinical routine in patients undergoing diagnostic cardiac catheterization or percutaneous cardiac interventional procedure performed with 7F-sheath.	50 consecutive patients Indication – diagnostic or invasive cardiac catheterization from common femoral artery	Primary endpoints - Primary efficacy endpoint – TTH (time to hemostasis) with or without oozing but without any additional compression, manual or mechanical. - Safety endpoint - incidence of adverse events - minor and major vascular complications and all other adverse events – serious and non-serious during duration of the study (i.e., until patient discharge).

PM-06723 FemoSeal™ VCS SSCP – Patients

1.4.3 Safety

FemoSeal™ has a shorter time to stop the bleeding at the procedure site compared to manual pressure. FemoSeal™ and manual pressure have similar risks. The clinical benefits of using FemoSeal™ outweigh the risks. Your doctor has considered the risks and decided that this device is appropriate for your procedure. See Section 2.1.4 for patient population information. If you have an existing health condition, speak with your doctor.

There were three post-market clinical studies for FemoSeal™. The 3 studies included a 2006 study with 100 patients, a 2008 study with 50 patients, and a 2009 study with 25 patients. These studies confirmed the safety and effectiveness of FemoSeal™. These studies also confirmed that FemoSeal™ performs as it is intended.

Terumo Medical Corporation continues to collect information on the device after its released in market. This is to demonstrate long-term safety and performance of the device. Terumo is conducting a study in Europe in 2022. This study will gather additional information for safety and performance of the FemoSeal™ device. In addition, a survey of physicians will be performed to gather information on the guidewire used with the device unit. The survey will gather information from physicians on what type of guidewire they use during the procedure.

1.5 Possible Therapeutic Alternatives

Manual pressure is the standard method to stop bleeding. It is the use of external hand pressure, by a member of the hospital staff, to the hole in the artery.

When considering other treatments, contact your doctor who can consider your medical history.

2.0 BIBLIOGRAPHY

Clinical Background Publications

1. Iannaccone M, Saint-Hilary G, Menardi D, et al. Network meta-analysis of studies comparing closure devices for femoral access after percutaneous coronary intervention. *J Cardiovasc Med* 2018; 19:586-596
2. Cox T, Lincourt A, Blair L, et al. Systematic Review of Randomized Controlled Trials Comparing Manual Compression to Vascular Closure Devices for Diagnostic and Therapeutic Arterial Procedures. *Surg Technol Int.* 2015; 27:32-44.
3. Tavriss DR, Wang Y, Jacobs S, et al. Bleeding and Vascular Complications at the Femoral Access Site Following Percutaneous Coronary Intervention (PCI): An Evaluation of Hemostasis Strategies. *J Invasive Cardiol.* 2012 Jul;24(7):328-34.
4. Jiang J, Zou J, Ma H, et al. Network Meta-analysis of Randomized Trials on the Safety of Vascular Closure Devices for Femoral Arterial Puncture Site Hemostasis. *Sci Rep.* 2015; 5: doi: 10.1038/srep13761

Publications Providing Safety and Performance Data on FemoSeal™ VCS

5. Alonzo, A., S. Rigattieri, F Giovannelli, et al. Transfemoral approach with systematic use of FemoSeal closure device compared to transradial approach in primary angioplasty. *Catheter Cardiovasc Interv* 2016; 87: 849-854.
6. Brenner, C., J. Margreitter, A. Gratl, et al. Femoral access site closure without prior femoral angiography: A retrospective analysis. *Wien Klin Wochenschr* 2018; 130: 197-203
7. Chivot, C., H. Deramond, R. Bouzerar and T. Yzet. Safety and Efficacy of Femoral Artery Closure with the FemoSeal Device After Cerebral Thrombectomy Using an 8 French Sheath. *Eur J Vasc Endovasc Surg* 2018; 55: 30-7344.
8. Choo, H. J., H. W. Jeong, J. Y. Park, et al. Ultrasonographic features of vascular closure devices: initial and 6-month follow-up results. *Ultrasonography* 2014; 33:283-290
9. De Poli, F., P. Leddet, P. Couppie, et al. FemoSeal Evaluation Registry (FER). Prospective study of femoral arterial closure with a mechanical system on 100 patients who underwent angioplasty procedures. *Ann Cardiol Angeiol (Paris)* 2014; 63:339-344.
10. Gabrielli, R., M. S. Rosati, M. Millarelli, et al. Device Use for Femoral Artery Closure by Different Techniques. *Ann Vasc Surg* 2018; 51:18-24

11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEvices Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019: doi:10.1177/1129729819854593

3.0 REVISION HISTORY

SSCP revision number	Date issued	Change description	Revision validated by the Notified Body
1	1-20-2023	Initial Issue, upon MDR CE certification	X Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)

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11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEVICES Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019: doi:10.1177/1129729819854593

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1. Iannaccone M, Saint-Hilary G, Menardi D, et al. Network meta-analysis of studies comparing closure devices for femoral access after percutaneous coronary intervention. J Cardiovasc Med 2018; 19:586-596
2. Cox T, Lincourt A, Blair L, et al. Systematic Review of Randomized Controlled Trials Comparing Manual Compression to Vascular Closure Devices for Diagnostic and Therapeutic Arterial Procedures. Surg Technol Int. 2015; 27:32-44.
3. Tavriss DR, Wang Y, Jacobs S, et al. Bleeding and Vascular Complications at the Femoral Access Site Following Percutaneous Coronary Intervention (PCI): An Evaluation of Hemostasis Strategies. J Invasive Cardiol. 2012 Jul;24(7):328-34.
4. Jiang J, Zou J, Ma H, et al. Network Meta-analysis of Randomized Trials on the Safety of Vascular Closure Devices for Femoral Arterial Puncture Site Hemostasis. Sci Rep. 2015; 5: doi: 10.1038/srep13761
5. Iannaccone M, Saint-Hilary G, Menardi D, et al. Network meta-analysis of studies comparing closure devices for femoral access after percutaneous coronary intervention. J Cardiovasc Med 2018; 19:586-596

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5. Alonzo, A., S. Rigattieri, F Giovannelli, et al. Transfemoral approach with systematic use of FemoSeal closure device compared to transradial approach in primary angioplasty. Catheter Cardiovasc Interv 2016; 87: 849-854.
6. Brenner, C., J. Margreitter, A. Gratl, et al. Femoral access site closure without prior femoral angiography: A retrospective analysis. Wien Klin Wochenschr 2018; 130: 197-203
7. Chivot, C., H. Deramond, R. Bouzerar and T. Yzet. Safety and Efficacy of Femoral Artery Closure with the FemoSeal Device After Cerebral Thrombectomy Using an 8 French Sheath. Eur J Vasc Endovasc Surg 2018; 55: 30-7344.
8. Choo, H. J., H. W. Jeong, J. Y. Park, et al. Ultrasonographic features of vascular closure devices: initial and 6-month follow-up results. Ultrasonography 2014; 33:283-290
9. De Poli, F., P. Leddet, P. Couppie, et al. FemoSeal Evaluation Registry (FER). Prospective study of femoral arterial closure with a mechanical system on 100 patients who underwent angioplasty procedures. Ann Cardiol Angeiol (Paris) 2014; 63:339-344.

10. Gabrielli, R., M. S. Rosati, M. Millarelli, et al. Device Use for Femoral Artery Closure by Different Techniques. *Ann Vasc Surg* 2018; 51:18-24
11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEVICES Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019; doi:10.1177/1129729819854593
16. Alonzo, A., S. Rigattieri, F. Giovannelli, et al. Transfemoral approach with systematic use of FemoSeal closure device compared to transradial approach in primary angioplasty. *Catheter Cardiovasc Interv* 2016; 87: 849-854.
17. Brenner, C., J. Margreitter, A. Gratl, et al. Femoral access site closure without prior femoral angiography: A retrospective analysis. *Wien Klin Wochenschr* 2018; 130: 197-203
18. Chivot, C., H. Deramond, R. Bouzerar and T. Yzet. Safety and Efficacy of Femoral Artery Closure with the FemoSeal Device After Cerebral Thrombectomy Using an 8 French Sheath. *Eur J Vasc Endovasc Surg* 2018; 55: 30-7344.
19. Choo, H. J., H. W. Jeong, J. Y. Park, et al. Ultrasonographic features of vascular closure devices: initial and 6-month follow-up results. *Ultrasonography* 2014; 33:283-290
20. De Poli, F., P. Leddet, P. Couppie, et al. FemoSeal Evaluation Registry (FER). Prospective study of femoral arterial closure with a mechanical system on 100 patients who underwent angioplasty procedures. *Ann Cardiol Angeiol (Paris)* 2014; 63:339-344.

21. Gabrielli, R., M. S. Rosati, M. Millarelli, et al. Device Use for Femoral Artery Closure by Different Techniques. *Ann Vasc Surg* 2018; 51:18-24
22. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
23. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEVICES Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
24. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
25. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
26. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019: doi:10.1177/1129729819854593

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2.0 BIBLIOGRAFIJA

Publikacije o kliničkoj pozadini

1. Iannaccone M, Saint-Hilary G, Menardi D, et al. Network meta-analysis of studies comparing closure devices for femoral access after percutaneous coronary intervention. J Cardiovasc Med 2018; 19:586-596
2. Cox T, Lincourt A, Blair L, et al. Systematic Review of Randomized Controlled Trials Comparing Manual Compression to Vascular Closure Devices for Diagnostic and Therapeutic Arterial Procedures. Surg Technol Int. 2015; 27:32-44.
3. Tavriss DR, Wang Y, Jacobs S, et al. Bleeding and Vascular Complications at the Femoral Access Site Following Percutaneous Coronary Intervention (PCI): An Evaluation of Hemostasis Strategies. J Invasive Cardiol. 2012 Jul;24(7):328-34.
4. Jiang J, Zou J, Ma H, et al. Network Meta-analysis of Randomized Trials on the Safety of Vascular Closure Devices for Femoral Arterial Puncture Site Hemostasis. Sci Rep. 2015; 5: doi: 10.1038/srep13761

Publikacije koje pružaju podatke o sigurnosti i učinkovitosti za sustav FemoSeal™ VCS

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10. Gabrielli, R., M. S. Rosati, M. Millarelli, et al. Device Use for Femoral Artery Closure by Different Techniques. *Ann Vasc Surg* 2018; 51:18-24
11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEVICES Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019: doi:10.1177/1129729819854593

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1. Iannaccone M, Saint-Hilary G, Menardi D, et al. Network meta-analysis of studies comparing closure devices for femoral access after percutaneous coronary intervention. J Cardiovasc Med 2018; 19:586-596
2. Cox T, Lincourt A, Blair L, et al. Systematic Review of Randomized Controlled Trials Comparing Manual Compression to Vascular Closure Devices for Diagnostic and Therapeutic Arterial Procedures. Surg Technol Int. 2015; 27:32-44.
3. Tavaris DR, Wang Y, Jacobs S, et al. Bleeding and Vascular Complications at the Femoral Access Site Following Percutaneous Coronary Intervention (PCI): An Evaluation of Hemostasis Strategies. J Invasive Cardiol. 2012 Jul;24(7):328-34.
4. Jiang J, Zou J, Ma H, et al. Network Meta-analysis of Randomized Trials on the Safety of Vascular Closure Devices for Femoral Arterial Puncture Site Hemostasis. Sci Rep. 2015; 5: doi: 10.1038/srep13761

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5. Alonzo, A., S. Rigattieri, F Giovannelli, et al. Transfemoral approach with systematic use of FemoSeal closure device compared to transradial approach in primary angioplasty. Catheter Cardiovasc Interv 2016; 87: 849-854.
6. Brenner, C., J. Margreitter, A. Gratl, et al. Femoral access site closure without prior femoral angiography: A retrospective analysis. Wien Klin Wochenschr 2018; 130: 197-203
7. Chivot, C., H. Deramond, R. Bouzerar and T. Yzet. Safety and Efficacy of Femoral Artery Closure with the FemoSeal Device After Cerebral Thrombectomy Using an 8 French Sheath. Eur J Vasc Endovasc Surg 2018; 55: 30-7344.
8. Choo, H. J., H. W. Jeong, J. Y. Park, et al. Ultrasonographic features of vascular closure devices: initial and 6-month follow-up results. Ultrasonography 2014; 33:283-290
9. De Poli, F., P. Leddet, P. Couppie, et al. FemoSeal Evaluation Registry (FER). Prospective study of femoral arterial closure with a mechanical system on 100 patients who underwent angioplasty procedures. Ann Cardiol Angeiol (Paris) 2014; 63:339-344.
10. Gabrielli, R., M. S. Rosati, M. Millarelli, et al. Device Use for Femoral Artery Closure by Different Techniques. Ann Vasc Surg 2018; 51:18-24

11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. Cardiovasc Intervent Radiol 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEvices Used in everyday Practice (CLOSE-UP) study. EuroIntervention 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. JAMA 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. Cardiovascular and interventional radiology 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. J Vasc Access 2019: doi:10.1177/1129729819854593

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11. Heo, Y. J., H. W. Jeong, J. W. Baek, et al. Ultrasound Evaluation of Puncture Sites After Deployment of Two Different Types of Vascular Closure Devices: A Prospective Comparative Study. *Cardiovasc Intervent Radiol* 2018; 41:1654-1663.
12. Holm, N. R., B. Sindberg, M. Schou, et al. Randomised comparison of manual compression and FemoSeal vascular closure device for closure after femoral artery access coronary angiography: the CLOSure dEVICES Used in everyday Practice (CLOSE-UP) study. *EuroIntervention* 2014; 10:183-190
13. Schulz-Schüpke, S., S. Helde, S. Gewalt, et al. Comparison of vascular closure devices vs manual compression after femoral artery puncture: the ISAR-CLOSURE randomized clinical trial. *JAMA* 2014; 312:1981-1987.
14. Sirignano, P., W. Mansour, A. d'Adamo, et al. Early experience with a new concept of angioplasty nitinol-constrained balloon catheter (Chocolate™) in severely claudicant patients. *Cardiovascular and interventional radiology* 2018; 41:377-384
15. Tagliaferro, F. B., G. Orgera, L. Mascagni, et al. FemoSeal™ vascular closure device for antegrade common femoral artery access: Safety and technical notes. *J Vasc Access* 2019; doi:10.1177/1129729819854593

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